TECHNICAL CHARACTERISTICS
Hopkins® THERMOFOCUS® Touch-Free Thermometer #571500 Series: Thermofocus 01500 Mod.: Thermofocus 01500A3em
Non-contact infrared thermometer

Resolution: 0.1

Body temperature readings			
Measurement range: 93.2-108.5° F (34.0-42.5°C)			
Room temperature working range: 50-104.0°F (10-40°C) ⁽¹⁾			
Accuracy °F	from 93.2 to 96.7°F:	±0.5°F	
	from 96.8 to 102.2°F:	±0.4°F ⁽¹⁾	
	from 102.3 to 108.5°F:	±0.5°F	
Accuracy °C	from 34.0 to 35.9°C:	±0.3°C	
	from 36.0 to 39.0°C:	±0.2°C ⁽¹⁾	
	from 39.1 to 42.5°C:	±0.3°C	

Other readings			
Measurement range: 33.8-131°F (1.0-55.0°C)			
Room temperature working range: 50-104.0°F (10-40°C) ⁽²⁾			
Accuracy °F	from 33.8 to 67.9°F	±1.8°F	
	from 68.0 to 96.7°F	±0.5°F	
	from 96.8 to 102.2 °F	±0.4°F	
	from 102.3 to 108.5°F	±0.5°F	
	from 108.7 to 131 °F	±1.8°F	
Accuracy °C	from 1.0 to 19.9°C:	±1.0°C	
	from 20.0 to 35.9°C:	±0.3°C	
	from 36.0 to 39.0°C:	±0.2°C	
	from 39.1 to 42.5°C:	±0.3°C	
	from 42.5 to 55°C:	±1.0°C	

⁽¹⁾ The ASTM E1965-98-2009 standard requires ±0.4°F (±0.2°C) accuracy for infrared thermometers in the 98.6-102.2°F (37-39°C) temperature range while the ASTM E567-86 and E112-86 standards require ±0.2°F (±0.1°C) precision for mercury and electronic thermometers in that same temperature range.

⁽²⁾ With room temperature between 50 and 60.8°F (10-16°C), if you use the "face" button, and with a room temperature between 5.5 and 10°C (41-50°F) if you use the "home" button, accuracy and the operating range are not guaranteed and the display shows "Lo.5" and the value alternately.

Power supply:	4 AAA alkaline batteries - 1.5 V
Life of high-quality batteries:	up to 3 years or 10,000 readings (depending on use)
Dimensions:	m 165 x 40 x 22 (6.5 x 1.45 x 0.86 inches) - including cap
Weight:	gr. 92 (3.17 oz.) - batteries included

Distance from the subject: **calculated using an optical aiming system** (approximately 1.2 inches/3 cm).

Large, clearly visible display.

Atmospheric pressure range of operating conditions: from 700 hPa to 1,060 hPa.

Relative humidity range of operating conditions: from 15 % to 93 %, non-condensing.

Keep in a clean, dry place, preferably at a temperature ranging between 60.8 and 104°F (+16 and +40°C). Store at a temperature included between -14 and 140°F (-10 and +60°C) and in any case not lower than 0.4°F (-18°C) or higher than 158°F (+70°C).

Expected life: 10 years.

The Hopkins® THERMOFOCUS® Touch-Free Thermometer is a class II medical device (II as per 93/42/EEC Directive and subsequent amendments and integrations) tested in hospitals, private clinics and medical offices. Contact the manufacturer for clinical tests.

The Hopkins® THERMOFOCUS® Touch-Free Thermometer is compliant with the pertinent ASTM standard (E1965-98-2009) and CEI EN 60601-1, CEI EN 60601-1-11 standards. Internally powered equipment for continuous operation.

The aiming LEDs of The Hopkins® THERMOFOCUS® Touch-Free Thermometer emit low luminous radiation largely under the level foreseen from IEC 60825-1:1993+A1:1997+A2:2001 Standard for Class I lights (class I, wavelength: 610 nm; maximum power: Pmax<1mW).

	Follow instructions for use
	ATTENTION: read the warnings
	Device designed for continuous use
	Applied part: type BF
	Recyclable material
IP 22	Protected against insertion of objects >0.5 inches (12.5 mm) & against vertical dripping water

  Dispose of the device in compliance with applicable legislation regarding electrical equipment and batteries.

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EN

USER MANUAL

Read the instructions carefully before using the thermometer

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Distance from the subject: **calculated using an optical aiming system** (approximately 1.2 inches/3 cm).

Large, clearly visible display.

Atmospheric pressure range of operating conditions: from 700 hPa to 1,060 hPa.

Relative humidity range of operating conditions: from 15 % to 93 %, non-condensing.

Keep in a clean, dry place, preferably at a temperature ranging between 60.8 and 104°F (+16 and +40°C). Store at a temperature included between -14 and 140°F (-10 and +60°C) and in any case not lower than 0.4°F (-18°C) or higher than 158°F (+70°C).

The Hopkins® THERMOFOCUS® Touch-Free Thermometer is compliant with the pertinent ASTM standard (E1965-98-2009) and CEI EN 60601-1, CEI EN 60601-1-11 standards. Internally powered equipment for continuous operation.

The aiming LEDs of The Hopkins® THERMOFOCUS® Touch-Free Thermometer emit low luminous radiation largely under the level foreseen from IEC 60825-1:1993+A1:1997+A2:2001 Standard for Class I lights (class I, wavelength: 610 nm; maximum power: Pmax<1mW).

1-800-835-1995
 www.hopkinsmedicalproducts.com

6850 Southbelt • Caledonia, MI 49316
 Cod. 915A6/H- r.00 141124

INDICATION FOR USE

The Hopkins® THERMOFOCUS® Touch-Free Thermometer is an infrared thermometer intended for intermittent measurement of human body temperature in people of all ages.

1. INTRODUCTION

Dear Customer, thank you for buying **The Hopkins® THERMOFOCUS® Touch-Free Thermometer.**

The Hopkins® THERMOFOCUS® Touch-Free Thermometer is truly easy to use. It is capable of measuring a child or adult's temperature without ever coming into contact with the skin: just bring it close to the forehead, the thermometer itself will tell you the correct distance. The Hopkins® THERMOFOCUS® Touch-Free Thermometer does not need to be placed inside any part of the body, and can even be used while a patient is sleeping.

2. WARNINGS

Read carefully before using The Hopkins® THERMOFOCUS® Touch-Free Thermometer

2.1 Precautions

- Use the thermometer in a room that is draft-free, at a steady temperature between 50 and 104°F (10 and 40°C).
- If the thermometer was stored in another room (or even in a drawer, cabinet, etc.) run the manual calibration procedure (MQCS) before using it. Alternatively, without touching the unit, wait at least 5 minutes for the temperature to stabilize or until the countdown on the display has been completed (see par. #8).
- Do not take a temperature reading if the person is sitting in a draft or if the subject has:
 - Been walking, running or exercising;
 - Come from another room that was ventilated or at a different temperature than the room where the thermometer is used;
 - Had a cool sponge or compress on their forehead;
 - Been wearing a cap, hat or scarf;
 - Been exposed to agents that could alter forehead temperature, e.g. a shower, shampoo, hair-drier, sponging, etc. (Even touching the forehead may alter the temperature.)
- In all the above cases, wait a few minutes for the forehead temperature to stabilize.
- Changing the reading point will give different results. Therefore, always aim the two lights on the same spot, precisely **at the center of the forehead** (midway between the top of the nose and the hairline) and keep the thermometer **perpendicular** to the forehead. Do not take measurements on areas other than the center of the forehead, except in the case of sweating (see par. #4.2)
- The temperature reading is taken in the area with a radius of ¼ inch (2 cm) around the point of light. It is of major importance to make sure that this area does not include eyebrows, hair or clothing. If necessary, brush away any hair from the forehead, but remember that this must be done a few minutes beforehand or the temperature reading will be higher than the real body temperature.
- When taking a temperature reading, remember that if oils, make-up, or an oxygen mask are present, and in the case of elderly, the temperature detected may be lower than the real body temperature.

- The forehead temperature reading can be affected by profuse sweating, superficial wounds or head injury.
- Do not use the thermometer on a sweaty forehead**, since the temperature reading would be unreliable. Read par. #4.2.
- The gilt waveguide (see figure 1) is the most delicate part of the device. It is composed of a gold-plated concave mirror that must be kept clean, crystal clear, and intact. Any damage, dust or dirt will alter the temperature reading.
- Do not handle the thermometer, and in particular the tip (fig. 1), for longer than strictly necessary before taking the reading.
- Do not use the thermometer in direct contact with the ear or any other part of the body.
- Do not use the thermometer in direct contact with objects or liquids. Do not submerge the thermometer in water or other liquids and keep it away from sources of heat and out of direct sunlight. If water seeps into the thermometer, contact Hopkins Medical Products® immediately for technical service.
- To prevent dust or dirt infiltrations in the tip, always close the cap.
- Do not use the device on a subject making a call with a mobile or cordless telephone or in the presence of strong electromagnetic fields.
- Avoid knocking or dropping the thermometer, and do not use it if it is damaged or not functioning properly.

2.2. Attention

1. Failure to observe the above-mentioned precautions may lead to very low or very high temperature readings, which cannot be attributed to product malfunctions.

- The aiming lights meet the photo-biological safety requirements.
- The unit is a delicate measurement instrument and must not be used by small children. It is not a toy. Keep it out of the reach of children or persons with limited sensorimotor skills. Small parts can be ingested or inhaled.
- Using this thermometer may not substitute medical consultation. Tell your physician what type of thermometer you are using and where on the body the temperature reading was taken.
- If the patient's movement makes it difficult to take a correct reading, turn the aiming lights on before bringing the thermometer close to the forehead.
- When taking your own temperature, use a mirror; or, if you have a child over six years of age, you can teach him to use The Hopkins® THERMOFOCUS® Touch-Free Thermometer.
- Since it never comes into contact with the body, The Hopkins® THERMOFOCUS® Touch-Free Thermometer does not require any "disposable" protective covers.

3. HOW IT WORKS

The Hopkins® THERMOFOCUS® Touch-Free Thermometer detects the infrared radiation coming from the human body. The forehead is an ideal site for taking a temperature reading, as it is crossed by the temporal artery and in direct contact with the brain. Finally, the head is the first part of the body to change temperature as a fever rises and falls. With each temperature measurement, in just a few hundredths of a second your Hopkins® THERMOFOCUS® Touch-Free Thermometer takes a series of 125 readings. Its sophisticated microprocessor then amplifies and processes this information along with the room temperature and shows the correct body temperature on the display.

It is important to remember that there is not a single "normal" temperature valid for everyone: moreover, individual temperature varies according to the measurement's site and throughout the day, also in response to physical or mental effort (for example, a crying child). Moreover, the body temperature can be affected by the outside temperature and, depending on the type of reading taken, other factors may also come into play.

Due to heat dispersion from uncovered parts of the body, the actual temperature at the forehead is generally lower than that in covered zones. Therefore, the Hopkins® THERMOFOCUS® Touch-Free Thermometer software automatically applies a correction factor and thus the resulting value is comparable to that given by other more usual temperature reading sites.

The Hopkins® THERMOFOCUS® Touch-Free Thermometer temperature reading taken on the forehead of a healthy person can range between 95 and 99.5°F (35 and 37.5°C).

To correctly judge a fever, readings must be taken under the same conditions each time (time of day, environment, etc.)

4. HOW TO USE THE THERMOMETER

4.1 Body Temperature

- Open the protective cap by rotating it 90° (fig. 2).
- Press the "face"  button and hold it down. The thermometer will switch on and the aiming lights will turn on.
- Hold the thermometer **perpendicular to the middle of the forehead** and move it closer to or further from the forehead until a single point of light appears. If the thermometer is too far away from the forehead, you will see two blurred points of light (fig. 3). If the thermometer is too close, you will see two separate points of light (fig. 4). When you see a **single point of light** on the subject's forehead, the thermometer is at the correct distance for an accurate temperature reading to be taken (fig. 5).
- Release the button and hold the device steady until the lights flash.
- Read the value shown on the display.
- If necessary, you may take another reading immediately.
- Close the protective cap.

The display will show the last reading for a few seconds, after which it will show the ambient temperature for 4 hours before turning itself off automatically.

4.2 In Case Of Sweating:

Reading On The Eyelid

At times, the forehead can show signs of sweat, for example as a fever drops. This can result in a low reading. Drying the forehead is not enough. In this case, you can take the temperature reading, still without contact, by taking a close-up scan of the closed eyelid (fig. 6). The light is harmless so there is no need to worry if the person who's temperature you are taking should open their eyes while you are taking a reading.

Precision is not guaranteed, but such reading can be considered a valid approximation of one's body temperature. In adults, such reading is also indicated when one has oil or make-up on the forehead; moreover it is also valid for the elderly.

4.3 Other Readings

The Hopkins® THERMOFOCUS® Touch-Free Thermometer can also be used to read the temperature of objects, liquids, and other surfaces in the 33.8-131°F (1-55°C) temperature range. For example:

- The temperature of a baby's feeding bottle (fig. 7), food, bath, etc. In this case, remember to mix the liquid well before taking the reading (with hot liquids or foods, take the reading quickly to prevent condensation from forming on the lens and wait 30 min. before taking another reading);
- Room temperature: by pointing the thermometer against a wall or a piece of furniture (see also par. #6);
- Temperature of an extremity, a wound, inflammation, or scar.

How to Perform Other Readings

Proceed as you would for a forehead reading but use the "home"  button.

In all these cases, if the taken temperature is lower than 93.2°F (34°C) or higher than 104.4°F (40°C), the given value will be displayed in alternation respectively with "Lo.3" or "Hi.2". These signs just inform the user that the taken values are outside the usual range of use (measurement of body temperature).

5. MEMORY FUNCTION

The memory function lets you call up the last 9 temperature readings. To activate this function, press the "Mem"  button: the value of the last reading will be shown on the display, accompanied by the number 1 and the symbol RE. Pressing the button again calls up the second to last, third to last reading and so on, accompanied by the numbers 2, 3, etc.

6. AMBIENT TEMPERATURE

When the thermometer is in standby mode ("Peak Performance" - see par. 7), the room temperature and the  symbol will flash and remain displayed for 4 hours after the last reading is taken.

If the display is off, press the "Mem"  button once and wait for 20 seconds; the room temperature will appear on the display.

7. HOW TO CHANGE

The Settings

Your thermometer leaves the factory preset to Fahrenheit (°F), and your thermometer can be preset also with display off (off  "Energy Savings") or display on (on  "Peak Performance") when in standby mode (4 hours after the last reading).

If necessary, these settings can be modified as follows:

- While the thermometer is off or in standby mode, press and hold the "Mem"  button; after about 8 seconds the display will change, showing the settings in sequence: °C; °F; "on"  ; "off" .
- When the desired setting appears, release the button. Only one setting can be modified at a time.

8. ROOM TEMPERATURE

CALIBRATION

If the device is handled at length or if there is a significant temperature difference with respect to the room temperature, a countdown may appear on the display. At this point you have the following options:

- Automatic Quick Calibration System (AQCS):** without touching the thermometer, wait until the countdown has run its course and the automatic quick calibration is completed. The countdown will continue updating as long as differences in temperature are detected (for example, because the thermometer has been handled). You can take a reading before the countdown has run its course, although in this case accuracy cannot be guaranteed, particularly if the temperature is rising.











- Manual Quick Calibration System (MQCS):** the alternative is to promptly correct the thermometer temperature, adapting it to the real temperature of the room where the reading is to be taken. Proceed as follows (the room temperature must be in the 50-104°F (10-40°C) temperature range):

- Press the "face"  and "home"  buttons fig. 8) simultaneously: the symbol CAL will appear on the display.
- Within 10 seconds open the cap and focus the thermometer on an internal wall (not the inside of an external wall) or wardrobe or cabinet with uniform temperature and at a point approximately 30 to 60 inches (80 to 150 cm) from the floor. Press the "home"  button (fig. 9).
- Release the button: the lights will flash slowly and the display will show the room temperature. To ensure a reliable temperature reading, do not focus the thermometer on an outside wall, window, source of heating or cooling (radiator, air conditioner, lamp, computer, a surface in contact with the human body, etc.).
- The thermometer is ready to take a reading.

Manual quick calibration (MQCS) can also be performed without the countdown if, for example, you need to move between rooms at different temperatures.


Both AQCS and MQCS enable the thermometer to take sufficiently accurate readings. When "RE" appears on the display, indicating that AQCS or MQCS have been performed.

- Normal calibration:** to achieve optimal readings, or to perform clinical tests, let the thermometer stabilize on its own. This requires waiting 10-30 minutes, depending on the difference in temperature.

9. MEANING OF DISPLAY MESSAGES			
ITEM DISPLAYED	DESCRIPTION	PROBLEM	REMEDY
	During normal operation, the battery symbol appears on the display.	The batteries are low but many measurements may be taken still.	Obtain new batteries for changing them when the signal E.3 (below) appears.
	"E.1" and the battery symbol appear on the display, or the unit does not turn on at all.	The batteries are dead.	Remove the batteries immediately and replace them when necessary (see par. #11).
	"Hi.4" appears on the display.	The room temperature is too high (above 104°F/40°C).	Move to another, cooler site.
	The display shows "Lo.5" and the value alternately.	The room temperature is between 50 and 60.6°F (10-15.9°C)	The temperature reading can be taken but accuracy is not guaranteed.
	The display reads "Lo.5".	The room temperature is too low (under the device's operative range).	Move to another, warmer room.
	The display shows "Hi.2" and the value alternately.	Temperature of 104°F (40°C) or more. It indicates a very high fever.	Seek medical attention.
	The display reads "Hi.2".	The body temperature appears too high.	Make sure that the instructions have been followed.
	The display shows "Lo.3" and the value alternately.	The surface temp. exceeds the limit for unit operation (>131°F/55°C).	The temperature reading cannot be taken.
	The display reads "Lo.3".	If a body temperature is being taken, this means that the body temperature appears to be too low (<93.2°F/34°C).	Check that the guide wave is not dirty or damaged, and that the subject has not just come from a location with a lower temperature.
	The display reads "Lo.3".	The body temperature appears too low.	Make sure that the instructions have been followed.
	The display shows a countdown, in minutes and seconds.	The device is stabilizing at the ambient temperature (AQCS).	Wait for automatic calibration AQCS, without touching the device, or run the MQCS (par. 8).
	The room temperature displayed is accompanied by the symbol RE.	The thermometer has undergone an automatic or manual calibration (AQCS or MQCS).	Precision is guaranteed, but to ensure a perfect reading, wait until normal thermometer calibration has been completed (10/30 minutes).
	The display reads "E.6".	The temperature of the thermometer is changing too fast.	Wait for the calibration of the device.
	The display reads "E.8".	The thermometer was moved before the light started flashing, or the area is subject to strong electromagnetic fields.	Wait until the light starts flashing before moving the thermometer; make certain that there are no mobile or cordless telephones in the vicinity.

10. TROUBLESHOOTING

- If the display does not turn on:
 - The batteries are completely dead or incorrectly inserted; replace or reinsert them (par. #11).
- If two bright points of light are visible on the subject's forehead:
 - The distance from the device to the forehead is incorrect; move the thermometer closer or further away until you see a single bright point of light (fig. 5).
- If the point of light is hard to see:
 - There is too much light in the room; cast a shadow over the subject.
- If the waveguide (fig. 1) is damaged or water has seeped into the thermometer:
 - Contact Hopkins Medical Products® for technical service.
- If the thermometer temperature reading is too low:
 - Make certain that the conditions outlined in the warnings (par. #2) have been met.
 - Check that the waveguide (fig. 1) is not soiled or damaged; if it is, clean it as indicated in par. #12 or contact Hopkins Medical Products® for technical service..
 - Check that the thermometer is perpendicular to the forehead as indicated in figure 5.
- If the thermometer temperature reading is too high:
 - Make certain that the conditions outlined in the warnings have been met (par. #2).
- If the thermometer appears blocked, does not turn off/does not go on standby or the aiming LEDs remain on after the button is released:
 - Reset the thermometer by removing and reinserting the batteries.

11. REPLACING THE BATTERIES

- Use your fingernail as shown in figure 10, or a pen, and press inside the hole on the back side of the device. Or insert the edge of a small coin in the crevice and apply a gentle leverage.
- Remove the battery cover.
- Remove the old batteries and dispose of them per your local regulations.
- Insert 4 new AAA batteries, preferably alkaline. Be sure to position the batteries as indicated in their housing.
- Close the battery cover.
- After changing the batteries, let the thermometer stabilize for 20 minutes before taking a temperature reading or perform the MQCS.
- **Remove the batteries if you do not expect to use the thermometer for a long period of time.**

12. CLEANING

CLEANING THE WAVEGUIDE

The thermometer waveguide is very delicate. Therefore, when the thermometer is not being used, we recommend that you always keep the cap on. However, if you need to remove dust or dirt from the waveguide or sensor at its base, use a cotton swab that has been slightly dampened with alcohol. Remove all dirt and make certain that nothing accumulates at the bottom of the waveguide where the sensor is located. Do not use any other objects or liquids as the surface of the waveguide or sensor could easily be scratched or damaged. Never let any excess liquid penetrate into the waveguide and sensor.

CLEANING THE THERMOMETER BODY

Use a soft cloth dampened with soap and water. Re-wipe with a sodium hypochlorite disinfectant if desired.

DO NOT USE the thermometer for at least 30 minutes after cleaning.

WARRANTY

Hopkins Medical Products® guarantees this product against any lack of conformity for 24 months as of the date of purchase (indicated on the your receipt or other fiscal document). Follow all instructions for the technical service procedures. This warranty does not cover the batteries and any damage caused by defective or run-down batteries or damage to the casing due to carelessness or improper use. The warranty is also voided if:

- The product is tampered with, damaged or used improperly.

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