

Apron Options & Sizing Order Form

Contact Name:


Date:

Email:

Phone:

Ship to Address:

Wrap-Around, Vest/Skirt Combo, or Adjust-A-Fit Aprons



Sizing Your Wrap-Around, Vest/Skirt Combo, or Adjust-A-Fit Aprons

Chest: Distance around chest at underarm

Waist: Distance around stomach at widest point

Hips: Distance around hips at widest point

(Size Charts available online)



WABR

VS

AFVS

Available Models

All Models have Patch Pocket only
Skirts have additional Slit Pocket

WABR.....Wrap Around Back Relief with Velcro

VSVest Skirt Combo with Velcro

AFVSAdjust-A-Fit Skirt with Velcro

STEP 1 Choose Model (ex: Vest Skirt Combo =VS) _____

STEP 2
Measure and place sizes into chart

	Wrap Around WABR	Vest VS, AFVS	Skirt VS, AFVS
Chest			N/A
Waist			
Hips		N/A	
Length Recommended length: top of knee	From sternal notch to desired length	From sternal notch to 2" below waist	Waist to desired length

Attention: Aprons are made to order and are non-returnable/non-refundable.
Please be sure measurements are accurate.

STEP 3 Choose Protection Option TL, TLL = Lite Ply 15-20% Less Weight (\$ Additional charges apply.)

☐ TL (0.5mm LE) ☐ TLL (0.35mm LE)

STEP 4 Choose Apron Fabric and Color Please Enter 2-3-Digit Code from page 3

STEP 5 Choose Belt Size — Back Relief (WABR only) Size: ☐ Small, 28"-38" ☐ Medium, 36"-48" ☐ Large, 44"-58"

STEP 6 Patch Pocket (Model #45576) ☐ Yes ☐ No (Left chest placement only) (\$ Additional charges apply).

STEP 7 Choose Apron Embroidery — Model # 5086107 (\$ Additional charges apply).

Embroidery is left chest placement only. **Block** Style. Up to 20 characters per line, 2 lines Maximum.

Please note: Embroidery is case sensitive and will be created exactly as written on lines below.

Choose Embroidery Color ☐ Black ☐ White
See page 2 for placement diagram

STEP 8 Add Logo (Additional charges apply) ☐ No Logo
Logo is left chest placement only. See page 2 for placement diagram.

- ☐ Customer Logo — Please send JPEG artwork with this form to orders@coneinstruments.com
☐ Vendor Stock Design — View options at www.shieldingintl.com/embroidery
☐ Medtronic Logo (Model #CCN17515)

STEP 9 Add Optional Thyroid Collar and Tether (\$ Additional charges apply.)

- Collar Options**
☐ Small Comfort Fit Thyroid Collar #5086193
☐ Large Comfort Fit Thyroid Collar #5086194

Add Tether*
☐

Choose Thyroid Collar Fabric and Color

Please Enter 2-3-Digit Code from page 3

*Collar tethered to apron on right side by a ring for easy removal


Apron Options & Sizing Order Form

Contact Name:	Date:
Email:	Phone:
Ship to Address:	

Apron Personalization (optional) (\$ Additional charges apply.)

STEP 10

Please indicate "L" for Logo, and/or "E" for Embroidery Placement Below:




Logo Placement = L
Embroidery Placement = E

Thyroid Collar Personalization (optional) (\$ Additional charges apply.)

STEP 11

Please indicate "L" for Logo, and/or "E" for Embroidery Placement Below:




Logo Placement = L
Embroidery Placement = E

Choose Thyroid Collar Embroidery* ☐ Black ☐ White

*16 characters maximum on Right and Left Sides Only
For small Comfort Fit Thyroid collar (#5086103) cannot be embroidered on the left strap. Right and front side only.

FABRICS & COLORS

Diamond Ripstop



Nylon



Designer Fabric

